

Radionuclides

Carbon-14 Handling Precautions

¹⁴C 5730 y β- 0.156 No γ E 0.156

Physical Data

Maximum Beta Energy: $0.156 \text{ MeV} (100\%)^{(1)}$ Maximum Range of Beta in Air: $22 \text{ cm} (8.6 \text{ in.})^{(2)}$

Occupational Limits(3)

Annual Limit on Intake: 2 mCi (74 MBq). Derived Air Concentration: 1 x 10^{-6} μ Ci/ml (37 kBq/m³).

Dosimetry

Millicurie (37 MBq) quantities of ¹⁴C do not present a significant external exposure hazard because the low-energy betas emitted barely penetrate the outer dead layer of skin. ¹⁴C-labeled compound uptake may be assumed to be uniformly distributed throughout all organs and tissues in the body⁽⁴⁾. Most ¹⁴C-labeled compounds are rapidly metabolized and the radionuclide is exhaled as ¹⁴CO₂. Some compounds and their metabolites are eliminated via the urine. Biological half lives vary from a few minutes to 40 days⁽⁴⁾.

PerkinElmer Life Sciences has developed the following suggestions for handling Carbon-14 after years of experience working with this low-energy beta emitter.

General Handling Precautions for Carbon-14

- 1. Designate area for handling ¹⁴C and clearly label all containers.
- 2. Prohibit eating, drinking, smoking and mouth pipetting in room where ¹⁴C is handled.
- 3. Use transfer pipets, spill trays and absorbent coverings to confine contamination.
- 4. Handle potentially volatile compounds in ventilated enclosures.
- 5. If enhanced containment is necessary, handle volatile compounds in closed systems vented through suitable traps.

- 6. Sample exhausted effluent and room air by drawing a known volume through a membrane filter followed by an impinger containing dilute NaOH.
- 7. Wear disposable lab coats, wrist guards and gloves for secondary protection.
- 8. Select gloves appropriate for chemicals handled.
- 9. Maintain contamination and exposure control by regularly monitoring and promptly decontaminating gloves and surfaces.
- 10. Use pancake or end-window Geiger-Mueller detectors or liquid scintillation counter to detect ¹⁴C.
- 11. Submit periodic urine and breath samples (as appropriate) for bioassay to determine uptake by personnel.



- 12. Isolate waste in sealed, clearly labeled containers and dispose according to approved guidelines.
- 13. Establish air concentration, surface contamination and bioassay action levels below regulatory limits. Investigate and correct any conditions that may cause these levels to be exceeded.
- 14. On completing an operation, secure all ¹⁴C; remove and dispose of protective clothing and coverings; monitor and decontaminate self and surfaces; wash hands and monitor them again.

Some 14 C-labeled compounds may penetrate gloves and skin. Handle these compounds remotely, wear two pairs of gloves and change the outer layer frequently. Special caution should be observed when handling 14 C-labeled halogenated acids. These compounds can be incorporated in the skin and deliver local dose commitments in the order of 10-100 rad per μ Ci (3-30 Gy per MBq) deposited.

References

- Kocher, David C., Radioactive Decay Data Tables, Springfield: National Technical Information Service, 1981 DOE/TIC-11026.
- 2. Kaplan, Irving, Nuclear Physics, New York: Addison-Wesley, 1964.
- U.S. Nuclear Regulatory Commission. 10CFR 20 Appendix B Standards for Protection Against Radiation, 1994.
- ICRP Publication 30, Part 3, Limits for Intakes of Radionuclides by Workers. Pergamon Press, Oxford, 1981.

This document contains general information designed to provide a basic understanding of radiation safety. While we believe the information to be accurate, regulatory requirements may change and information contained herein is not tailored to individual needs. A radiation protection specialist should be consulted for specific applications.



Worldwide Headquarters: PerkinElmer Life Sciences, 549 Albany Street, Boston, MA 02118-2512 USA (800) 551-2121

Technical Support: in Europe: perkinelmer.europe@perkinelmer.com in US and Rest of the World: techsupport@perkinelmer.com