


Laboratory Waste disposal Flow chart




Tiger stripe tagged Yellow caddy

Yellow caddy

Tiger stripe tagged Yellow caddy

Autoclaves which reach 134°C are available in room 2.42 in corridor 3 on the 2nd floor and room 1.24 in corridor 4 on the 1st floor Jenner Wing.





Place the box in a Yellow lidded caddy beneath the above sign or a red tagged yellow caddy Place the box or bag in a Yellow lidded caddy beneath the above sign or a purple tagged yellow caddy

Place the waste bag in a Yellow lidded caddy beneath the above sign or a Tiger stripe tagged yellow caddy

Autoclaves which reach 134°C are available in room 2.42 in corridor 3 on the 2nd floor and room 1.24 in corridor 4 on the 1st floor Jenner Wing.